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
Chief Operating Officer
Giampiero Belcredi

Organismo Notificato n. 0476
Notified Body nr. 0476

SANISIM SOLUTION 6 % H₂O₂

DATA SHEET
Rev. 3 del 12/09/2017

DATA SHEET SANISIM SOLUTION

1	Applicant's name and address	
	Society	Dimensione Service s.a.s. : Via III Trav Ludovico D'Angiò, 22 70032 Bitonto (BA) - ITALIA
	Telephone	: +39 080-3752362
	Fax	: +39 080-2143172
	e-mail address	: info@dimensioneservice.it
2	STATEMENT OF PRODUCT (see annex certificate)	
	Active principle	
	Chemical name	: SANISIM SOLUTION 6% H ₂ O ₂
	CAS Nr.	: 7722-84-1
	CEE Nr.	: 008-003-00-9
	Molecular formula	: H ₂ O ₂
	Other component	
	Chemical name	: Silver nitrate
	CAS Nr.	: 7761-88-8
	CEE Nr.	: -
	Molecular formula	: AgNO ₃
	Quantitative formula	
	hydrogen peroxide	: 6%
	Silver nitrate	: 40,8mg/l
	distilled water	: 100% c.s.p.

3	PHYSICS APPEARANCE	
	SANISIM SOLUTION It is presented in liquid form	

4	PRODUCT USE AND JUSTIFICATION USE
----------	--

The system allows SANITISING all surfaces including those that are generally inaccessible or hidden in hospitals or other hazardous environments such as crowded. The tests and performed experiments show a great disinfectant activity on microorganisms. SANISIM is based on the synergy of two principles: the hydrogen peroxide more silver cations: (H₂O₂ + Ag⁺).

The mechanisms of action are:

- 1) The oxidizing action of hydrogen peroxide which generates:
 - Oxidation of lipid membrane structures.
 - An alteration of ribosomes and nucleic acids.
- 2) Cathodic Action of Ag⁺, implies:
 - A inversion of polarity of membrane with a consequent alteration.
 - An inhibition of synthesis protein.

These two combined actions ensure the elimination of microorganisms. The system uses a spreading and spraying technology.

Ionization:

The sprinkler system produces a quality aerosol mixed with droplets of size about 8µm, equipped with electrostatic charge to facilitate the nucleation process.

SANISIM SOLUTION 6 % H₂O₂

DATA SHEET
Rev. 3 del 12/09/2017

Nucleation:

Nucleation is a process by means of which small droplets adhere both to atmospheric particles on the surfaces, creating an interface between the droplets and particles on one side, and between the substance in the gaseous phase and the liquid disinfectant from other side.

The SANISIM method ensures a higher degree of disinfection, because the disinfectant fully is present on all the surfaces to be treated. The method is the only one that provides to an activity of broad-spectrum disinfection.

5	PROPERTY 'PHYSICAL AND CHEMICAL	
	Physical and chemical properties	Value
	Denomination	SANISIM SOLUTION
	Physical state	Liquid
	Smell	Slightly acidic
	Colour	Colorless
	Density at 20 ° C (kg/l)	0,9646
	pH a 20°C	4,0 ± 0,5
	Fusion point (°C)	-55
	Boiling point (°C)	~ 110
	Temperature of autodecomposition (°C)	111
	Inflammability piont	No data available
	Explosive Limit	No data available
	Auto ignition point or spontaneous ignition	No data available
	Water solubility	Total
	Viscosity at 20°C (cP)	0,98
	Vapor pressure at 20°C (mmHg)	18,30
	Saturation of vapor in air at 25°C (mg/m ³)	787

6	Saturation of vapor in air at	
	Product is distributed in bottles in PE from 1L to be used in SANISIM device. To transport the product it may be contained in tanks of 3L, 5L and 20L.	
	Distribution spraying	: Sanitization Full 6 ml/m ³
	Distribution spraying	: Sanitization preventive 3 ml/m ³

7	DURATION OF PRODUCT (see the stability annex tests)
	Product is an oxidizing agent and responsive, stable if you follow the recommended conditions of storage and handling. The commercial product is stabilized to prevent decompositions for contamination. SANISIM SOLUTION stabilized for 1 year from date of manufacture. The products are stabilized to prevent the contaminants to keep him losing his business.

8	PACKAGING AND SALE
	HDPE plastic bottles with inner cap and cap sealed 1L

9	STATEMENT OF PRODUCTION FACTORY
	The currently laboratory that produces small amounts of SANISIM SOLUTION è:
	Raro Srl – Industria Detergenti Professionali
	Via 1° Maggio, 14 – 75100 MATERA (MT)
	Tel. 0835 383370 – Fax 0835 383473
	MILANO – Tel. 02/96460832 – Fax 02/96460048
	VAT: CE IT 00116260779
	e.mail: info@rarosrl.com
	Internet site: http://www.rarosrl.com/

SANISIM SOLUTION 6 % H₂O₂

DATA SHEET
Rev. 3 del 12/09/2017

	PRODOTTO
	N ° of registration in Register of medical devices class IIb CE n° 0546
10	CERTIFICATE ANALYSIS OF ACTIVE PRINCIPLES AS-QUANTITY (See attached)
11	SAMPLES OF DIFFERENT PRODUCTION BATCH NR. or LOT WILL BE TAKEN FROM THE LABORATORY AND MADE AN ASSESSMENT OF CHEMICAL PHYSICS TO CONFIRM THE VALIDITY 'OF PRODUCT (to do periodically to qualify SANISIM solution manufacturer)
12	ANALYTICAL METHODS TO DETERMINING ACTIVE INGREDIENT
	To verify the contents and wealth of hydrogen peroxide proceeds to value the active component of the product:
	1. SANISIM solution it is diluted to a concentration of 0.5-0.6%
	2. In 10 ml, It is added 1 g of KI diluted in 100ml of water and 20ml of H ₂ SO ₄ (1V / 2V of water) shaking all
	3. Iodine that will develop is titrated with sodium thiosulfate 10N until the disappearance of the brown color. You can add starch to improve the loss of color.
	<i>1ml sodium thiosulfate is equivalent to 1.7001 mg of hydrogen peroxide</i>
	Note: To accelerate the reaction it can add a 2% ammonium molybdate
	valuation method of silver nitrate method:
	The silver nitrate is incorporated in the same formula. The product is very pure and chlorides with absence of any other impurities, for which it is sufficient that the solution in distilled water is total, leaving no residue and a completely transparent and colorless. To analyze the other active ingredient of the product, on: atomic absorption spectrophotometry with silver lamp.
13	CHEMICAL AND PHYSICS INFORMATION ON STABILITY
	EXPLOSIVE PROPERTY
	The dissolutions of H ₂ O ₂ up to 70% are not considered explosive. However, explosion may take place when the hydrogen peroxide is mixed with organic compounds to form a single phase, an emulsion or a suspension. Determining factors include: <ul style="list-style-type: none"> - Concentrations with a mixture of H₂O₂, H₂O and organic matter; - type of organic matter - presence an initiator - temperature of the mixture - Presence of a decomposition catalyst in the mixture. H ₂ O ₂ It is highly explosive when it is placed in contact with strong oxidizing agents or reducers.
	COMBURENT PROPERTY
	It is not an explosive product under normal conditions. It can be broken down by heat, radiation or impurities, which cause explosive gases or intervenes in oxidation reactions with heat dispersion and risk of combustion.
	CORROSIVE PROPERTY
	When using as a disinfectant, the product is always used in a watery solution. The result of corrosion of materials of evidence suggests no or only few corrosive effects of watery solution of product of common materials, such as: steel, aluminum, wood and polymers.
	PHYSICAL CHEMISTRY INCOMPATIBILITY WITH OTHER PRODUCTS
	Incompatibility: High temperatures, metals, oxidants and reducers, acids, alkalis, halogens, electromagnetic radiation, organic solvents, etc. In the case of decomposition, it breaks down into oxygen and water according to the following equation::

SANISIM SOLUTION 6 % H₂O₂

DATA SHEET
Rev. 3 del 12/09/2017

	<p style="text-align: center;">$2\text{H}_2\text{O}_2 \rightarrow 2\text{H}_2\text{O} + \text{O}_2 + 98 \text{ kJ per mole}$</p> <p>In an alkaline the rate of decomposition rapidly increases with the increase of pH. The H₂O₂ absorbs a broad spectrum of radiation, for this light can cause a decomposition photochemistry. The heating of the hydrogen peroxide to external agents is important. On the basis of physical-chemical criteria, the rate of decomposition both at the surface as in solution can duplicate itself every 10 ° C increase in temperature. In case of a contamination of the reaction is accelerated. The most dangerous contaminants are the salts of iron, copper, vanadium, tungsten, molybdenum, silver or platinum. This phenomenon is known with the name of homogeneous decomposition. Small amounts of contaminants, may originate a rapid decomposition. The effect of pH on the contaminated peroxide decomposition rate is very high.</p> <p>al theories on the mechanism of decomposition reaction of H₂O₂.</p> <p>metal ions, it comes to chain reactions in which said ions oxidize and reduce subsequently, a small amount for this reason can cause decomposition of a large amount of peroxide.</p> <p>When it is placed in contact with the hydrogen peroxide can be produced a solid material insoluble rapid decomposition. the hydrogen peroxide decomposes on contact with any surface at room temperature. The rate of decomposition varies greatly according to the material and the state of said surface. The oxides and hydroxides of heavy metals and the noble metals are very active decomposition catalysts. The activity is greater the greater the specific surface area.</p>
--	---

14	EFFICIENCY (See annex evidence)	<p>Thanks to its great destructive power of germs and its magnificent residual behavior, you can spray SANISIM solution as a disinfectant to wide range disinfectant for all type of surfaces and environments:</p> <p>It is effective against:</p> <ul style="list-style-type: none"> - Sterilization: bacteria Gram + e Gram -. - Virulicidal : wide range - Fungicide : wide range. - Sporicide : bacterial spores and fungi. - Protozoa: amoeba, giardia, etc. - Oocysts of coccidia: Eimeria spp, cryptosporidium. - Alghicida: algae. <p>The solution SANISIM activity follows the rules of the effectiveness of disinfectants European Committee for Standardization: EN 1275, EN 1276, EN 1040 and EN 1650.</p> <p>Amply exceeds the requirements of these requests the European Standards on chemical antiseptics and disinfectants.</p>
	UNI EN 1040:2006	Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of basic bactericidal activity of chemical disinfectants and antiseptics - Test method and requirements (phase 1).
		<i>Pseudomonas aeruginosa; Staphylococcus aureus, Escherichia coli; Enterococcus sp.</i>
	NF T 72-281	Determination of bactericidal, fungicidal, sporicidal activity for aerial surface disinfection processes
	UNI EN 13697:2001	Chemical disinfectants and antiseptics - Quantitative Test for non-porous surfaces for the evaluation of bactericidal and / or fungicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas - Test method and requirements without mechanical action (phase 2 / step 2)
		"Chemical disinfectants and antiseptics - Preservation of test

SANISIM SOLUTION 6 % H₂O₂

DATA SHEET
Rev. 3 del 12/09/2017

	UNI EN 14353:2013	organisms used for determining of bactericidal (including Legionella), mycobactericidal, sporicidal, fungicidal and virucidal (including bacteriophages)"
	UNI EN 13727:2014	Chemical disinfectants and antiseptics - Quantitative suspension test for evaluation of bactericidal activity in medical area - Test method and requirements (phase 2, step 1)
	UNI EN 14561:2006	Chemical disinfectants and antiseptics - Quantitative carrier test for evaluation of germs of bactericidal activity for instruments used in medical area - Test method and requirements (phase 2, step 2)
	UNI EN 14562:2006	Chemical disinfectants and antiseptics - Quantitative carrier test for the evaluation of fungicidal or yeasticidal activity for instruments used in the medical area - Test method and requirements (phase 2, step 2)
	UNI EN 14563:2009	Chemical disinfectants and antiseptics - Quantitative carrier test for the evaluation of mycobactericidal or tuberculocidal activity of chemical disinfectants used for instruments in the medical area - Test method and requirements (phase 2, step 2)
	UNI EN 14348:2005	Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of mycobactericidal activity of chemical disinfectants in the medical area including instrument disinfectants - Test methods and requirements (phase 2, step 1)
	UNI EN 13704:2005	Chemical disinfectants - Quantitative suspension test for the evaluation of sporicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas - Test method and requirements (phase 2, step 1)
	The effectiveness of the disinfectant in front of bacteria and fungi is also demonstrated in the information of the evidence made by the head of the laboratory Dr. SERINO TECNOLAB place at premises of our company. The full feature set of efficacy and final report reports for all bacteria considered.	

15	TOXICITY PRODUCT (See annex evidence)
	<p>Toxicity tests made SANISIM solution of acute toxicity, eye irritation and skin irritation are in according to Directive 67/548 / EEC. The tests made with a disinfectant dose more than twice that recommended by the manufacturer, have given the following results:</p> <ul style="list-style-type: none"> - Study Eye irritation: SANISIM solution does not produce; - Study of skin inflammation: SANISIM solution does not produce.

THE RECOMMENDED DOSES, SANISIM COLUTION IS TOTALLY HARMLESS TO MEN AND ANIMALS AND TO ALL AREAS WHERE TO APPLY.

16 SAFETY DATA SHEET (See attached)

17 CLASSIFICATION AND LABELING OF PRODUCT IN THE LAW (See attached)
Copies of labels (See attached)

18 INDICATIONS AND QUANTITATIVE OF ACTIVE EXPRESSED IN g per 100g (See annex Certificate of Analysis)

19 INSTRUCTIONS AND METHOD USE (See Label)

20 WARNINGS AND SAFETY PHRASES

SANISIM SOLUTION 6 % H₂O₂

DATA SHEET
Rev. 3 del 12/09/2017

SAFETY STANDARDS IN HANDLING, MIXING AND APPLICATION

- For handling it is recommended to place the product in accessible places, hydrants of water to proceed as quickly as possible, the dilution of the product in the event of an accident.
- In the place of handling you must have good ventilation. You must have showers and washrooms. II Staff should be trained on the hazards of the product and the forms of embodiment in the case of an accident.

TREATMENT METHOD AND ELIMINATION OF PRODUCT IN EXCESS AND ITS CONTAINERS

The disposal of the product will be carried out by neutralization with water. The containers and packaging will not be used for other products.

SAFETY TIME FOR ENTRY IN AREAS TREATED

Contact time 15 minutes before entering into the environment with any required Personal Protective Equipment.

SPECIAL PRECAUTIONS IN STORAGE AND TRANSPORT:

The types of containers where one has to store the product will be of aluminum or stainless steel. The containers can be of PVC, PE, PTFE. The tanks will be closed for collection pens with anti leak channeling. The product is incompatible with metals and plastics. The product should be stored in cool and ventilated places, separated from incompatible products and high temperature sources

Transport information:

ONU Number: Not classified

Risk class: CORROSIVE

IMO/IMDG: Not classified

IMDG: Not classified

ADR/RID: Not classified

Packaging group: Not classified

N ° identification of the danger: Not classified

PANELS ORANGE NUMBERED: Not classified

PRECAUTIONS IN CASE OF FIRE. NATURE OF GAS IN BURNING

It is not a product it will flammable fuel and its solutions, under normal conditions, are not explosive. It can decompose by heat, radiation or metal impurities which give rise to explosive vapors or intervene in oxidation reactions with loss of heat or fire risk.

Means of extinction: Use a thin layer of water to absorb gases and films to cool equipment, containers, etc.

Inadequate methods of extinction: none known

Special instructions against fire: it will use water spray to cool containers exposed to the fire still not caught fire, and to dilute the missing products. You must use appropriate clothing and respiratory protection due to product vapors. Cool with water containers and stores as they may explode for gas formation. Working backs to the wind. Extraordinary danger to fire: Risk of corrosive vapors of the product and dispersion of oxygen (oxidizing).

Special protective clothing in case of fire: self-contained breathing clothes and adequate

OTHER TOOL OF SAFETY

TOOL TO BE IMPLEMENTED IN ACCIDENTAL LOSS:

Personal precautions: keep staff not protected far from infected area with wind behind. Avoid contact with losses and use appropriate protective clothing.

Precautions for products: use a film of water to absorb gases of fumes. Prevent that the product Log in to the water mains tubes, tanks, closed rooms. Dilute with plenty of water.

Implementation methods: signal the area and follow current legislation on the environment.

SANISIM SOLUTION 6 % H₂O₂DATA SHEET
Rev. 3 del 12/09/2017**FIRST AIDS**

Danger nature: corrosive action on live tissues. Risk of suffocation and permanent lung injury.

basic information about the first aid: remove person concerned from danger zone. Wash with plenty of water, keep still and fix it. Get medical attention.

Effects of exposure: burning, irritation of mucous membranes, coughing.

Actuations for Eyes: Wash with plenty of water for 15 minutes at least and go to the eye doctor urgently.

Clothing: remove them immediately, showering.

Skin: wash the affected area with plenty of water and consult a doctor.

Ingestion: If conscious, give water and keep it repaired. Do not induce vomiting. keep immobile.

Inhalation: Move the affected person from the contaminated area, in the open air, and practice artificial respiration if necessary.

Need for medical care: Always urgent.

Are needed in workplaces: Showers and security sinks, eye and respiratory protection masks gloves, suits and boots resistant to acids

NOTE: With sending of product is attached safety data sheet

RESIDUAL

This information is not necessary because product SANISIM solution is applied as a dry fog and spray on all surfaces to be treated. It is a totally biodegradable product, does not leave any kind of residue.

good Company

If you have any questions, please contact

Dr. Venanzio Bielli

President - E2S3 Srl

+39.329-4527646

venanzio.bielli@e2s3.it

vb_e2s3

<https://www.linkedin.com/in/venanzio-bielli-336b9830>

vbielli-e2s3

+39.329.4527646

Italian Headquarters
E2S3 Srl
Via Volta n° 3
20028 Legnano (Milan)
Italy

China branch - 中国分公司

E2S3 co, Ltd
大连市中山区友好路 227 号悦泰城里 412,
412 Yuetai City Lane, No.227 Youhao
Road, Distretto di Zhongshan, Dalian

Poland Branch Office
Przemyslaw Kaminski
East Europe Strategies Director &
Board Member
Paczków
Telephon: +48-774339445
Mobile: +48-606200267
przemyslaw.kaminski@e2s3.eu

孙妍丽 – Sun Yanli
亚洲项目发展经理和董事会成员
Manager for Asia projects' development
& Board Member

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